□□□ LEADERGUN®
Made in Italy

User Manual

User Manual

rev. 0.3 of 12-05-2017

ATTENTION:

Please, read carefully this manual before using Leadergun

Congratulations!

Thank you for purchasing

the most innovative injector for Mesotherapy

Completely designed and manufactured in Italy


Mesotech Aps

Juelsmindevej 57, 7120 Vejle, Denmark

www.mesotech.dk

Description

Main features

LEADERGUN® is one of the most advanced fully automated injection system.

LEADERGUN® has been developed thanks to Mesotech long experience in mesotherapy and has been designed to fix all the problems connected to portable devices.

LEADERGUN® is completely microprocessor-controlled and thanks to its color display (D) lets you to constantly view the programming values set.

LEADERGUN® is the only portable device which, thank to the use of a microprocessor controlled by stepper motor, can make accurate and precise injections without medicine leakage.

LEADERGUN® is able to work with 1.0ml, 2.0 ml, 2.5ml, 5.0ml and 10.0ml standard mesotherapy syringes and 4.0mm, 6.0mm and 13.0mm needles without any external adapter kit. The needle penetration depth is completely automatic and can be set, via a microprocessor controlled by stepper motor, between 1.0mm and 13.0mm.

LEADERGUN® has a Cover (C) used to protect the syringe from accidental operation during treatment.

LEADERGUN® is powered by an universal AC / DC adapter that provide power to the battery for recharging.

LEADERGUN® is equipped with a rechargeable high capacity Li-Ion battery that allows to operate without AC / DC adapter for several hours.

Package contents: N° 1 Charger Brand Globe Tek, Inc. mod. GTM91099-6015-3.0-T2, N° 1 Power Cord 100cm long, N° 5 Skin Stabilizer, N° 1 Instructions Manual, Warranty conditions.

Technical specifications

AC/DC adapter Power Box Medline EXM805118:

| | |
|--|----------------|
| Input voltage: | 100v – 240v ac |
| Frequency: | 50 / 60 Hz |
| Output voltage: | 12V dc |
| Max output current: | 5A dc |
| Power Plug Europe, UK, USA, New Zeland, Australia. | |

Li– Ion Battery:

| | |
|------------------|------------|
| Nominal Voltage: | 7.4V |
| Capacity: | 2200mA / h |

Storage conditions:

| | |
|----------------------------|-------------------|
| Minimum temperature: | -10°C |
| Maximum temperature: | 75°C |
| Minimum Relative Humidity: | 30% no-condensing |
| Maximum relative humidity: | 95% no-condensing |

Low atmospheric pressure: 700hPa
Maximum atmospheric pressure: 1060hPa

Operation:

Minimum temperature: 10° C
Maximum temperature: 40° C
Minimum relative Humidity: 30% non-condensing
Maximum relative humidity: 75% non-condensing
Low atmospheric pressure: 700hPa
Maximum atmospheric pressure: 1060hPa

Classification and Symbols:

Equipment not suitable for use in the presence of flammable anesthetics with air, oxygen or nitrous oxide

Device suitable for continuous operation

Device not protected against the penetration of liquids

Attention, check the documentation attached

Equipmentb subject to the provisions of the UE Directive (WEEE)

Device with Type B Applied Part

Device of Classe II b

Continue current

AC

Power Jack

Button (trigger) to initiate the dispen the product (G)

Manufactured by

IPX0


Maintenance and Technical Documentation:

In case of failure and / or request for technical specifications, please contact:

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Analysis test

THE DEVICE MUST BE USED IN ACCORDANCE WITH THE FOLLOWS EMC INFORMATIONS

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Intended use

LEADERGUN "mesotherapy injector" is an electronic portable device for fully automatic injecting systems, it allows accurate and accurate injections without waste of medicine.

Environment: The Device is mainly used in the medical field in outpatient environments.

Operator type: The device can be used only by doctors possessing the necessary know how, about the use of DM in the field of treatments.

Clinical Purpose of the Device: Portable Electronic Device for Fully Automatic Injection Systems.

Restrictions on the use of the device: Transport and Storage: Minimum Temperature: -10° C, Maximum Temperature: 75° C, Minimum Relative Humidity: 30% Non Condensing, Maximum Relative Humidity: 95% No Condensation, Atmospheric Pressure Minimum: 700hPa, Maximum Atmospheric Pressure: 1060hPa

Operating: Minimum Temperature: 10° C, Maximum Temperature: 40° C, Minimum Relative Humidity: 30% Non Condensing, Maximum Relative Humidity: 75% No Condensation, Atmospheric Pressure Minimum: 700hPa, Maximum Atmospheric Pressure: 1060hPa.

Warnings - Possible side effects

No undesirable effects on the use of the device are observed when used in accordance with the manufacturer's instructions for use and maintenance.

Contraindications

Limits of use: The Device must be used in accordance with the EMC information referred to in Section 02. Radio equipment and mobile phones, if used near the Device, may affect the operation of the device. Use of accessories and cables other than those specified, with the exception of those sold by the device manufacturer as spare parts for internal components, may cause an increase in emissions or a decrease in the immunity of the device.

The Device should not be used near other Devices if it is not, you must be careful that the Device will continue to operate in its normal configuration.

Potential Hazard Situations:

environmental: if used outside the parameters indicated in this user manual.

Clinical Patient: according to the judgment of the physician using the device

Regarding device conditions: the device must be clean; syringes, needles must be sterile and disposable; the skin stabilizer (accessory supplied with DM) should be single use.

Of interaction between different devices or other devices: the Device interacts with the human body through the contact established by the "stabilizer for the skin" (accessory supplied with DM) with the area to be treated. In addition, the device operates with other devices such as sterile and disposable sterile syringes and needles that inoculate the suitable drug for treatment.

Handling of the product

Precautions and Warnings: LEADERGUN must be used in accordance with the EMC information contained in Annex 1 referred to in the last pages of this manual.

Radio equipment and cellular phones, if used near the LEADERGUN, may affect the operation of the device.

Use of accessories and cables other than those specified, with the exception of those sold by the LEADERGUN manufacturer as spare parts for internal components, may cause an increase in emissions or a decrease in the immunity of the appliance.

LEADERGUN should not be used near other Appliances otherwise you should be careful that LEADERGUN will continue to function in its normal configuration.

Cleaning the LEADERGUN:

WARNING: Maintenance is recommended before and after each use. Carefully clean the parts to be used at each treatment. Before and after each use, clean the body of the LEADERGUN and especially the end of the appliance and the sliding part with a cloth dampened with aseptic product or cotton soaked with 90 ° alcohol.

Cleaning, Disinfection and Sterilization of Patient Contact Parties:

For more hygiene, sterile and single use syringes and needles should be used. Skin stabilizers (DM accessory) are disposable.

LEADERGUN Li-Ion Battery Management:

WARNING: Always charge the battery after a long period of inoperability.

For short periods of inoperability, it is advisable to leave the Charger Connected to the LEADERGUN Power Jack.

LEADERGUN penetration capacity:

In order to guarantee sufficient penetration capacity of the needle, it is advisable not to use LEADERGUN with the tip upwards. During use, the appliance must be held horizontally or with the tip down.

Device storage Mode:

In the original packaging (Carrying case).

Predictable improper use:

It is not foreseeable for improper use by the user if he or she strictly adheres to the manufacturer's instructions for use and maintenance of the DM.

Connecting and Disconnecting Leading Parts of LEADERGUN

Syringe Placement:

- Prepare the syringe with the product to be injected.
- Raise Cover (C).
- Turn on LEADERGUN.
- From the Home Home, by Joystick (J), select the desired work program (Glow, Lipo, Stretch, PRP or Botulinum) and press Enter.
- LEADERGUN will upload the selected work program.
- The Cute Stabilizer (A) will automatically move in such a way that it respects the programmed depth of penetration.
- Next, the Carrellino (B) will also move automatically and position it in such a way as to easily accommodate the syringe to be used.
- The Display (D) will display the PLEASE WAIT and then the Syringe Housing Improvement screen where the Joystick (J) up and down keys can refine the syringe housing.
- Place the needle in the Syringe Base Attachment (E).
- Next, place the syringe plunger in the Punching Punch Socket (F) compartment in the seat corresponding to the ml capacity of the syringe used.
- Get down the Cover (C).
- WARNING: Do not force the Push Tiller Cartridge (B) manually, but always use the Syringe Housing Enhancement Function.
- To enable this function, press the center button of the Joystick (J) for about 3s with the cursor positioned on SETUP.

Positioning of the Pocket Cute Stabilizer Port:

The LEADERGUN is equipped with a Metal Stabilizer Door Pad located at the front. It is used to accommodate the Cute Stabilizer (A) and define, by means of a stepper motor, the depth of the needle.

At power-up, once the work program is loaded, it will be positioned according to the preset depth data, while turning off will go back to its zero position.

Positioning of the Skin Stabilizer:

Once the Skin Stabilizer Push Point is positioned according to the preset depth data, you can accommodate the Skin Stabilizer (A) that serves to display the inoculation point of the product.

Connecting to the Power Supply

- Insert the power cord or plug of the Charger into an electrical outlet.
- Connect the Charger Output Cable to the Power Jack (K) located at the bottom of the device handle.
- The battery charging symbol will appear on the display.
- Wait about 3:00 hours for full charge or the Display shows battery charge.
- Battery charge is automatically handled by the microprocessor so an overload will not cause damage to the battery.
- The battery can not be switched on while charging the battery.

Removable parts and accessories used:

All accessories supplied with the LEADERGUN are referred to in the "Package Contents" section.

The use of other unqualified accessories by the manufacturer diminishes the safety of the appliance.

Their use demands responsibility to the user.

Functional features

How to recharge LEADERGUN battery:

Plug the AC/DC adapter's power cord into an electrical socket compliant with the specifications mentioned above.

Connect the AC / DC adapter's output cable in the Plug (J) located at the base of the device.

The display will show the battery charge icon.

Wait for about 3:00 hours to fully charge or the Display to give an indication of battery charged.

The battery charge process is automatically handled by the microprocessor so an overcharging will not cause damage to it.

While charging the battery you can turn the device on and possibly operate.

How to turn on the LEADERGUN:

Press for few seconds the middle Joystick (J) button.

The display lights on will shows LEADERGUN logo as shown in fig.1

Fig. 1. Turning on

Immediately after the Display (D) shows the main Home page as shown in fig.2

Fig. 2. Home

How to turn off the LEADERGUN:

Press for 3sec. the middle joystick (J) button with the cursor positioned on PW OFF. The LEADERGUN will turn off and if still plugged to the AC/DC adapter, a battery charge icon will be shown on the Display (D)

How to use the LEADERGUN:

You can select the following modes with the joystick (J), from the main Home page:

- **GLOW:** The LEADERGUN will load the MESOGLOW program in which working mode is preset by default as follows:
 - Syringe: 5.0ml
 - Needle: 4.0mm
 - Penetration depth: 2.0mm
 - Drop: micro
 - Penetration time: 0.2s
 - Injection: single (spot)
- or with the SETUP function you can set the following parameter
 - **Syringe: 1.0 ml, 2.0 ml, 2.5ml, 5.0ml, 10.0ml**

- Needle: 4.0mm, 6.0mm, 13mm
- Penetration depth: from 1.0mm to 13.0mm
- Drop:

| Syringe | Pico | Micro | Macro |
|---------|---------|---------|---------|
| 1.0ml | 0.007ml | 0.013ml | 0.025ml |
| 2.0ml | 0.013ml | 0.025ml | 0.05ml |
| 2.5ml | 0.013ml | 0.025ml | 0.05ml |
| 5.0ml | 0.025ml | 0.05ml | 0.10ml |
| 10.0ml | 0.05ml | 0.10ml | 0.15ml |

+/-20%

- Needle's time into the skin injection after:
 2.0ms, 3.0ms, 4.0ms, 5.0ms, 6.0ms, 7.0ms, 8.0ms, 9.0ms, 10.0ms
- Injection: single or repetitive
 In case of repetitive injection you can set the number of injections per minute. In details:

| Syringe | Pico | Micro | Macro |
|---------|------------|-----------|----------|
| 1.0ml | 150-120-80 | 116-82-55 | 82-51-34 |
| 2.0ml | 141-109-73 | 107-73-49 | 72-44-29 |
| 2.5ml | 141-109-73 | 107-73-49 | 72-44-29 |
| 5.0ml | 120-86-57 | 82-52-34 | 51-29-19 |
| 10.0ml | 91-59-39 | 55-32-21 | 40-22-14 |

+/-20%

- LIPO: The LEADERGUN will load the MESOLIPO program in which working mode is preset by default as follows:
- Syringe: 10.0ml
- Needle: 13.0mm
- Penetration depth: 9.0mm
- Drop: macro
- Needle's time into the skin injection after: 0.2s
- Injection: single (spot)
 or with the SETUP function you can set the following parameters:
- Syringe : 1.0 ml, 2.0 ml, 2.5ml, 5.0ml, 10.0ml
- Needle: 4.0mm, 13.0mm
- Penetration depth from: 1.0mm to 13.0mm

- Drop:

| Syringe | Pico | Micro | Macro |
|---------|---------|---------|---------|
| 1.0ml | 0.007ml | 0.013ml | 0.025ml |
| 2.0ml | 0.013ml | 0.025ml | 0.05ml |
| 2.5ml | 0.013ml | 0.025ml | 0.05ml |
| 5.0ml | 0.025ml | 0.05ml | 0.10ml |
| 10.0ml | 0.05ml | 0.10ml | 0.15ml |

+/-20%

- Needle's time into the skin injection after:
 2.0ms, 3.0ms, 4.0ms, 5.0ms, 6.0ms, 7.0ms, 8.0ms, 9.0ms, 10.0ms
- Injection: single or repetitive
 In case of repetitive injection you can set the number of injections per minute. In details:

| Syringe | Pico | Micro | Macro |
|---------|------------|-----------|----------|
| 1.0ml | 150-120-80 | 116-82-55 | 82-51-34 |
| 2.0ml | 141-109-73 | 107-73-49 | 72-44-29 |
| 2.5ml | 141-109-73 | 107-73-49 | 72-44-29 |
| 5.0ml | 120-86-57 | 82-52-34 | 51-29-19 |
| 10.0ml | 91-59-39 | 55-32-21 | 40-22-14 |

+/-20%

- STRETCH: The LEADERGUN will load the STRETCHMARKS program in which working mode is preset by default as follows:
 - Syringe: 10.0ml
 - Needle: 13.0mm
 - Penetration depth: 12.0mm
 - Flux: 1.0ml / min
 - Needle's time into the skin injection after: 0.4s
 - Injection: continuous
 or with the SETUP function you can set the following parameters:
 - Syringe : 1.0 ml, 2.0 ml, 2.5ml, 5.0ml, 10.0ml
 - Needle: 4.0mm, 13.0mm
 - Penetration depth from: 1.0mm to 13.0mm

- Injection: single or repetitive

In case of repetitive injection you can set the number of injections per minute. In details:

| Syringe | Pico | Micro | Macro |
|---------|------------|-----------|----------|
| 1.0ml | 150-120-80 | 116-82-55 | 82-51-34 |
| 2.0ml | 141-109-73 | 107-73-49 | 72-44-29 |
| 2.5ml | 141-109-73 | 107-73-49 | 72-44-29 |
| 5.0ml | 120-86-57 | 82-52-34 | 51-29-19 |
| 10.0ml | 91-59-39 | 55-32-21 | 40-22-14 |

+/-20%

- TOXIN The LEADERGUN will load the TOXIN program in which working mode is preset by default as follows:

- Syringe:1.0ml
- Needle: 13.0mm
- Penetration depth: 2.0mm
- Drop: micro
- Needle's time into the skin injection after: 0.4s
- Injection: single (spot)

or with the SETUP function you can set the following parameter:

- **Syringe : 1.0 ml, 2.0 ml, 2.5ml, 5.0ml, 10.0ml**
- Needle: 4.0mm, 6.0mm; 13.0mm
- Penetration depth from: 1.0mm to 13.0mm
- Drop:

| Syringe | Pico | Micro | Macro |
|---------|---------|---------|---------|
| 1.0ml | 0.007ml | 0.013ml | 0.025ml |
| 2.0ml | 0.013ml | 0.025ml | 0.05ml |
| 2.5ml | 0.013ml | 0.025ml | 0.05ml |
| 5.0ml | 0.025ml | 0.05ml | 0.10ml |
| 10.0ml | 0.05ml | 0.10ml | 0.15ml |

+/-20%

- Needle's time into the skin injection after:

2.0ms, 3.0ms, 4.0ms, 5.0ms, 6.0ms, 7.0ms, 8.0ms, 9.0ms, 10.0ms

| Syringe | 1.0ml / min | 5.0ml / min | 10ml / min |
|---------|-------------|-------------|------------|
| 1.0ml | OK | OK | N. A. |
| 2.0ml | OK | OK | OK |
| 2.5ml | OK | OK | OK |
| 5.0ml | OK | OK | OK |
| 10.0ml | OK | OK | N. A. |

- Needle's time into the skin injection after:
 2.0ms, 3.0ms, 4.0ms, 5.0ms, 6.0ms, 7.0ms, 8.0ms, 9.0ms, 10.0ms
 - Injection: continuous
 - LATEST: The LEADERGUN will return to the mode before the shutdown.
 - MR: The LEADERGUN will load from the internal memory a mode previously stored by the operator.
 - PRP: The LEADERGUN will load the PRP program in which working mode is preset by default as follows:
 - Syringe: 5.0ml
 - Needle: 13.0mm
 - Penetration depth: 2.0mm
 - Drop: micro
 - Needle's time into the skin injection after: 0.4s
 - Injection: single (spot)
- or with the SETUP function you can set the following parameter:
- **Syringe : 1.0 ml, 2.0 ml, 2.5ml, 5.0ml, 10.0ml**
 - Needle: 4.0mm, 6.0mm; 13.0mm
 - Penetration depth from: 1.0mm to 13.0mm
 - Drop:

| Syringe | Pico | Micro | Macro |
|---------|---------|---------|---------|
| 1.0ml | 0.007ml | 0.013ml | 0.025ml |
| 2.0ml | 0.013ml | 0.025ml | 0.05ml |
| 2.5ml | 0.013ml | 0.025ml | 0.05ml |
| 5.0ml | 0.025ml | 0.05ml | 0.10ml |
| 10.0ml | 0.05ml | 0.10ml | 0.15ml |

+/-20%

- Needle's time into the skin injection after:
 2.0ms, 3.0ms, 4.0ms, 5.0ms, 6.0ms, 7.0ms, 8.0ms, 9.0ms, 10.0ms

- Injection: single or repetitive

In case of repetitive injection you can set the number of injections per minute. In details:

| Syringe | Pico | Micro | Macro |
|---------|------------|-----------|----------|
| 1.0ml | 150-120-80 | 116-82-55 | 82-51-34 |
| 2.0ml | 141-109-73 | 107-73-49 | 72-44-29 |
| 2.5ml | 141-109-73 | 107-73-49 | 72-44-29 |
| 5.0ml | 120-86-57 | 82-52-34 | 51-29-19 |
| 10.0ml | 91-59-39 | 55-32-21 | 40-22-14 |

+/-20%

How to access the MESOGLOW program

From the main Home, select, with the joystick (J), the GLOW icon and press Enter.

The LEADERGUN will load the MESOGLOW program in which the icons are marked by a small green triangle.

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Fig. 3. GLOW– SETUP - SPOT

Once the program has been loaded, the Trolley (B) will move automatically and will position in such a way to easily accommodate the syringe loaded with 5.0ml.

The Display (D) will show first a “PLEASE WAIT” message and then the fig.4 screen in which it's possible to refine the syringe accommodation using UP and DOWN joystick button.

Fig. 4. Refining syringe accommodation

After about 5 sec., or pressing the middle Joystick (J) button, the Display (D) will show fig.3 screen once again.

In SETUP mode the Trigger (G) will be inactive while the Joystick (J) will be active in order to ensure maximum safety against accidental injections.

To activate the Trigger (G) and then operate according to preset parameters, you must switch from SETUP mode to OPERATE mode.

To do that, position the cursor on SETUP and press Enter. The Display (D) will show the following screen:

Fig. 5. GLOW – OPERATE - SPOT

Except for its middle button, the Joystick (J) will be inactive, while the Trigger (G) will be active. At each press of the Trigger (G) the LEADERGUN will make a single injection (SPOT) according to preset parameters.

How to accesss the MESOLIPO program

From the main Home, select, with the joystick (J), the LIPO icon and press Enter. The LEADERGUN will load the MESOLIPO program, in which the icons are marked by a small red triangle. The Display (D) will show the screen in fig.6

Fig. 6. LIPO - SETUP - SPOT

Once the program has been loaded, the Trolley (B) will move automatically and will position in such a way to easily accommodate the syringe loaded with 10ml. The Display (D) will show a “PLEASE WAIT” message and than the fig.7 screen in which it's possible to refine the syringe accommodation using UP and DOWN joystick button.

Fig. 7. Refining syringe accomodation

After about 5 sec., or pressing the middle Joystick (J) button, the Display (D) will show fig.6 screen once again.

In SETUP mode the Trigger (G) will be inactive while the Joystick (J) will be active in order to ensure maximum safety against accidental injections.

To activate the Trigger (G) and than operate according to preset parameters, you must switch from SETUP mode to OPERATE mode.

To do that, position the cursor on SETUP and press Enter.

The Display (D) will show the following screen:

Fig. 8. LIPO – OPERATE - SPOT

Needle's time into the skin injection after

Except for its middle button, the Joystick (J) will be inactive, while the Trigger (G) will be active. At each press of the Trigger (G) the LEADERGUN will make a single injection (SPOT) according to preset parameters.

How to access the STRETCHMARKS program

From the main Home, select, with the joystick (J), the STRETCH icon and press Enter.

The LEADERGUN will load the STRETCHMARKS program, in which the icons are marked by a small blue triangle.

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Fig. 9. STRETCH - SETUP

Once the program has been loaded, the Trolley (B) will move automatically and will position in such a way to easily accommodate the syringe loaded with 10ml.

The Display (D) will show a "PLEASE WAIT" message and then the fig.10 screen in which it's possible to refine the syringe accommodation using UP and DOWN joystick button.

Fig. 10. Refining syringe accommodation

After about 5 sec., or pressing the middle Joystick (J) button, the Display (D) will show fig.9 screen once again.

In SETUP mode the Trigger (G) will be inactive while the Joystick (J) will be active in order to ensure maximum safety against accidental injections.

To activate the Trigger (G) and then operate according to preset parameters, you must switch from SETUP mode to OPERATE mode.

To do that, position the cursor on SETUP and press Enter.

The Display (D) will show the following screen:

Fig. 11. STRETCH - OPERATE

Except for its middle button, the Joystick (J) will be inactive, while the Trigger (G) will be active. At each press of the Trigger (G) the LEADERGUN will make a single injection (SPOT) according to preset parameters.

Except for its middle button, the Joystick (J) will be inactive, while the Trigger (G) will be active. At each press of the Trigger (G) the LEADERGUN will make a single injection (SPOT) according to preset parameters.

How to access the TOXIN program

From the main Home, select, with the joystick (J), the TOXIN icon and press Enter.

The LEADERGUN will load the TOXIN program, in which the icons are marked by a small blue triangle.

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Fig. 9. TOXIN - SETUP

Once the program has been loaded, the Trolley (B) will move automatically and will position in such a way to easily accommodate the syringe loaded with 10ml.

The Display (D) will show a "PLEASE WAIT" message and then the fig.10 screen in which it's possible to refine the syringe accommodation using UP and DOWN joystick button.

Fig. 10. Refining syringe accommodation

After about 5 sec., or pressing the middle Joystick (J) button, the Display (D) will show fig.9 screen once again.

In SETUP mode the Trigger (G) will be inactive while the Joystick (J) will be active in order to ensure maximum safety against accidental injections.

To activate the Trigger (G) and then operate according to preset parameters, you must switch from SETUP mode to OPERATE mode.

To do that, position the cursor on SETUP and press Enter.

The Display (D) will show the following screen:

Fig. 11. TOXIN - OPERATE

Except for its middle button, the Joystick (J) will be inactive, while the Trigger (G) will be active. At each press of the Trigger (G) the LEADERGUN will make a single injection (SPOT) according to preset parameters.

Except for its middle button, the Joystick (J) will be inactive, while the Trigger (G) will be active. At each press of the Trigger (G) the LEADERGUN will make a single injection (SPOT) according to preset parameters.

How to access the PRP program

From the main Home, select, with the joystick (J), the PRP icon and press Enter.

The LEADERGUN will load the STRETCHMARKS program, in which the icons are marked by a small blue triangle.

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Fig. 9. PRP- SETUP

Once the program has been loaded, the Trolley (B) will move automatically and will position in such a way to easily accommodate the syringe loaded with 10ml.

The Display (D) will show a "PLEASE WAIT" message and then the fig.10 screen in which it's possible to refine the syringe accommodation using UP and DOWN joystick button.

Fig. 10. Refining syringe accommodation

After about 5 sec., or pressing the middle Joystick (J) button, the Display (D) will show fig.9 screen once again.

In SETUP mode the Trigger (G) will be inactive while the Joystick (J) will be active in order to ensure maximum safety against accidental injections.

To activate the Trigger (G) and then operate according to preset parameters, you must switch from SETUP mode to OPERATE mode.

To do that, position the cursor on SETUP and press Enter.

The Display (D) will show the following screen:

Fig. 11. PRP- OPERATE

Except for its middle button, the Joystick (J) will be inactive, while the Trigger (G) will be active. At each press of the Trigger (G) the LEADERGUN will make a single injection (SPOT) according to preset parameters.

How to change preset parameters

In any programs is possible to change preset parameters.

In SETUP mode select with the Joystick (J) the icon with parameter to be changed and press Enter.

The Display (D) will show new parameters by mean of the Left and Right Joystick (J) buttons.

Select the new parameter and press Enter to confirm.

Select another icon, if necessary, to repeat the operation or switch to OPERATE mode to operate with the new parameters.

How to set Display contrast

While changing parameters is possible to change Display (D) contrast with Up and Down Joystick (J) buttons.

How to work in Repetitive Injection mode

It's possible to work in repetitive (burst) injection mode only in MESOGLOW and MESOLIPO programs.

To do that, select SPOT icon in SETUP mode and press Enter.

The LEADERGUN will switch from SPOT to BURST mode and the display will show fig 12 screen if you are in MESOGLOW program or fig.13 if you are in MESOLIPO program:

Fig. 12. GLOW – SETUP – BURST

Fig. 13. LIPO – SETUP – BURST

In both cases is possible to change the number of injections per minute.

To change this setting select the middle icon with the Joystick (J) and press Enter.

The Display (D) will show new parameters by mean of the Left and Right Joystick (J) buttons.

Select the new parameter and press ENTER to confirm.

The screen with the new parameter will be shown on the Display(D).

Select another icon, if necessary, to repeat the operation or switch to OPERATE mode to work with the new parameters.

To do that, position the cursor on SETUP and press Enter

The Display (D) will show fig 14 screen if you are in MESOGLOW program or fig.15 if you are in MESOLIPO program:

Fig. 14. GLOW - OPERATE - BURST

Fig. 15. LIPO - OPERATE - BURST

Except for its middle button, the Joystick (J) will be inactive, while the Trigger (G) will be active. At each press of the Trigger (G) the LEADERGUN will make a cycle of repetitive injections (BURST) according to parameters set.

To return to single injection mode (SPOT), press Enter to switch to SETUP mode.

Select BURST function with the Joystick (J) and press Enter.

The LEADERGUN will return to single injection mode (SPOT).

How to store changed parameters in memory

In MESOGLOW, MESOLIPO and STRETCHMARKS is possible to store in memory the new configuration once you have changed the parameters.

Position the cursor on SETUP with the Joystick (J), select SAVE function and press Enter.

Select with Left and Right Joystick (J) buttons the memory location where you want to store the new configuration.

The Display (D) will show fig.16 screen if the selected memory location is still occupied otherwise will show fig.17 screen if memory location is empty.

Fig. 16. Save occupied

Fig. 17. Save empty

Press Enter to store in memory and return to the new working configuration.
To exit from SAVE function without storing the new configuration in memory, select EXIT (fig.18) with Left Joystick (J) button, then press Enter

Fig. 18. Memoria Exit

How to return to main HOME

In SETUP mode, position the cursor on BACK and press Enter.
The Display (D) will show the main Home screen (fig.2).

How to return to working program before shutdown

In the main Home select LATEST function with the Joystick (J) and press Enter.
The LEADERGUN will load the last program before shutdown and the Display (D) will show the parameters set.

How to recall from memory a saved work program

From main Home select MR function with the Joystick (J) and press Enter.
The Display (D) will show fig.19 screen concerning the first program stored in memory.

Fig. 19. MR 1

Select the desired work program stored using Left and Right Joystick (J) buttons and press Enter.

The LEADERGUN will load the program from the memory and the Display (D) will show the parameters set.

To exit from MR function without recalling a stored work program select EXIT (fig.20) using Left Joystick (J) button and press Enter


Fig. 20 MR Exit

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In this document, the Manufacturer guarantees and declares, under his own responsibility, that the device placed on the market satisfies the essential requirements of the Directive 93/42 / EEC and of the Legislative Decree n. 46 of 24/02/1997.

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Norme applicabili / *Applicable standards:*

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(*Legal representative of Mesotech S.r.l.*)

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Please, follow this link and complete the online form:

<http://www.salute.gov.it/DispoVigilancePortaleRapportoOperatoreWeb/>

Please also send the copying report for knowledge also to:

mesotechsrl@arubapec.it

needle
tech®

 needleflex®
NON-BRUIING EXTRA SMOOTH CANNULA

- Reduces the pain of the injection, without damaging the dermis and underlying tissues
- Needleflex penetrates gently thanks to a homogeneous lubrication present on its entire surface
- Wide range of sizes and lengths
- Suitable for every medical application: filler and mesococktail.
- Perfect exhibility and strength of the steel tube and the plastic "hub"
- Excellent proportion between quality and price
- Minimum tissue trauma due to the blunt tip and the flexibility of the cannula.
- Needleflex reduces the risk of bruising and inflammatory reactions.
- A box contains 100 pieces

| Cannula | Size | Length |
|---------|------|--------|
| Luer | 30G | 25mm |
| Luer | 27G | 38mm |
| Luer | 27G | 50mm |
| Luer | 25G | 38mm |
| Luer | 25G | 50mm |
| Luer | 22G | 50mm |
| Luer | 22G | 100mm |
| Luer | 18G | 50mm |
| Luer | 18G | 100mm |

- Needletech hypodermic needles are made with a careful selection of raw materials, maintaining a particular attention during the packaging phase, in order to avoid any defect of the tips of the needles.
- An exclusive lubrication process of the needle cannula allows a smooth and painless penetration into the skin.
- Accurate sharpening of needles minimizes the pain associated with the injection.
- Ultrathin diameters (0.3 and 0.20 mm) make this product excellent for mesotherapy, filler, botulinum and sclerosing treatments.
- 100 piece box.
- In sterile single-use EO blister packs.

| Cannula | Size | Length |
|---------|------|--------|
| Luer | 33G | 4mm |
| Luer | 32G | 8mm |
| Luer | 32G | 4mm |
| Luer | 31G | 13mm |
| Luer | 31G | 6mm |
| Luer | 31G | 4mm |
| Luer | 30G | 13mm |
| Luer | 27G | 13mm |
| Luer | 26G | 13mm |

Mesotech is specializing in providing aesthetic and skin care technologies for the beauty and medical institutes. Customer satisfaction and quality are the first priorities for our staff. We develop our Range in a continuous dialogue with our customers. Development of new products from original ideas to deep research. Original formulations where unique ingredients are selected to normalise and correct skin disorders and imbalances

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(*Rappresentante Legale della Mesotech S.r.l.*)